Práctico N° 1

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An efficient method for the extraction of DNA from vertebrate tissues

With the advent of PCR (Ref. 1), the extraction of DNA for genetic analysis has become one of the most time-consuming procedures in the laboratory, particularly when a large number of samples must be processed. Numerous DNA extraction procedures have been published. While some of these are undoubtedly rapid, most require either the use of hazardous chemicals or the purchase of exotic, often expensive reagents and equipment. We have adapted two existing methods^{2,3} to produce a simple and efficient DNA extraction procedure that uses common laboratory reagents. This method has been used to obtain pure, high molecular mass DNA from whole blood, skin and liver biopsies obtained from a variety of vertebrate species. Amplifiable DNA has also been obtained from plucked hair samples using this method.

Each sample (100 mg of tissue or 100 μ l of whole blood) was suspended in 1.5 ml microfuge tubes with 300 μ l of digestion buffer (100 mM NaCl, 50 mM Tris-HCl, 1% SDS, 50 mM EDTA, pH 8.0). Proteinase K (10 mg ml⁻¹) was added to a final concentration of 100 μ g ml⁻¹ and the samples incubated for 2 h at 50°C, then overnight at 37°C. After digestion was complete, an equal volume (300 μ l) of 5 M LiCl was added to each tube (5 M NaCl or 2.5 M NaClO₄ can be substituted, but the DNA quality and yield decreases slightly, see Table 1). The sample was mixed thoroughly by inversion for 1 min, then 600 μ l of chloroform added and the samples placed on a rotating wheel for 30 min. The sample was spun for 15 min at high speed in a benchtop microfuge and the supernatant carefully removed to a new microfuge tube. Exactly 2 volumes of room temperature absolute ethanol were added and the tube inverted several times until the DNA precipitated. The DNA was recovered by high speed centrifugation, the supernatant decanted, and the DNA pellet washed briefly with 70% ethanol. Excess ethanol was removed and the sample allowed to air dry for 10 min. The sample was then resuspended in 100–200 μ l of TE buffer (10 mM Tris-HCl, 1 mM EDTA, pH 7.5) and allowed to dissolve overnight at 4°C.

TABLE 1. Spectrophotometric comparison of DNA extraction procedures^a

	Phenol	NaCl	NaClO ₄	LiCl
No. of samples	4	4	4	34
Mean spectrophotometric absorbance ratio (260/280 nm)	1.79	1.61	1.75	1.79
Mean DNA yield (μ g) from 100 mg sample ^b	227	105	260	318

^aTissue samples were extracted as described in the text and the spectrophotometric absorbances determined at 260 and 280 nm.

^bDNA yields were extrapolated from the absorbance at 260 nm.